



# ComforTrac

comfortable • simple • effective

## Cervical Traction




The ComforTrac cervical home traction unit has been carefully engineered to provide maximum comfort, efficacy, and ease of use. Details that highlight this superior design include a push-of-a-button “quick release” pump, and a contouring memory foam headrest that cradles the head. Likewise, transporting the unit has been made effortless through the use of a lightweight durable carry-case.

# Cervical Traction

Our innovative design makes it possible for patients to go home with the same quality traction they receive during physical therapy.

**Fully Adjustable Angle of Incline** without the need to order extra parts




**Hand Pump**  
our patented hand pump makes it easy to increase or decrease traction force

**Quick Release Button**  
push to release the traction force

**Memory Pillow**  
conforms comfortably to the head

**Independently Adjustable Neck Wedges**  
provide a custom fit



## SPECIFICATIONS:

**Maximum force:** 50 lbs.

**Angles of incline:** 10, 15 and 20 degrees (with no additional parts or attachments)

**Hand pump:** Displays force in both pounds and kilograms

**Neck wedges:** Independently adjustable

**Headrest material:** Memory foam

**Testing:** 15,000 cycles which equals using 3 times daily for 15 years

## SUPPORTING RESEARCH AND REFERENCES:

Spine. 1996 Aug 15;21(16):1877-83. "Nonoperative management of herniated cervical intervertebral disc with radiculopathy." PMID: 8875719

Spine. 1985 Dec;10(10):867-71. "Cervical spine disorders. A comparison of three types of traction." PMID: 3914085

J Manipulative Physiol Ther. 2002 Mar-Apr;25(3):188-92. "Intermittent cervical traction for cervical radiculopathy caused by large-volume herniated disks." PMID: 11986581

Bangladesh Med Res Counc Bull. 2002 Aug;28(2):61-9. "Effects of cervical traction and exercise therapy in cervical spondylosis." PMID: 12825763



ComforTrac Inc.  
PO Box 1056  
McLean, VA 22101

703.891.0455

info@comfortrac.net  
[www.comfortrac.net](http://www.comfortrac.net)